

B'SYS GmbH

CHO K_v4.3 Cell Line

Specification Sheet

© B'SYS GmbH

TABLE OF CONTENTS

1. PRODUCT SHIPMENT	3
1.1. Product Format	3
1.2. Mycoplasma Certificate	3
1.3. Assays	3
2. VALIDATION OF CHO KV4.3 CELLS	4
2.1. Electrophysiology	4
2.2. Biophysical characterization	5
2.3. K_v4.3 channel Blocker Dapoxetine	5
2.4. Patch-clamp Success Rates	6
3. CELL CULTURE CONDITIONS	6
3.1. General	6
3.2. Recommended Complete Medium	6
3.3. Antibiotics	6
3.4. Thawing Cells	6
3.5. Splitting Cells	7
3.6. Freezing Cells	7
3.7. Stability of CHO K_v4.3 cells	7
4. KV4.3 SEQUENCE	7
4.1. Accession Number NM_002232	7
5. CONTACT INFORMATION	8
5.1. Contact Address for Technical Support and Ordering Information	8

1. PRODUCT SHIPMENT

1.1. Product Format

CHO cells stably transfected with recombinant K_v4.3 potassium channel:

- 0.75 mL aliquots of frozen cells at 1 E+06 cells/mL
- Cells are frozen in complete medium with 10% DMSO
- Cells are frozen at passage number, see vial

1.2. Mycoplasma Certificate

B'SYS periodically tests cells for presence of mycoplasma by means of highly sensitive PCR based assays. All delivered cells are free of mycoplasma.

1.3. Assays

CHO K_v4.3 cells were validated for manual patch-clamping and automated patch-clamping (Q-Patch, Sophion)

2. VALIDATION OF CHO KV4.3 CELLS

2.1. Electrophysiology

K_v4.3 currents were measured by means of the patch-clamp technique in the whole-cell configuration. The bath solution contained (in mM) NaCl 137, KCl 4, CaCl₂ 1.8, MgCl₂ 1, D-glucose 10, HEPES 10, pH (NaOH) 7.40. The pipette solution consisted of (in mM) KCl 130, MgCl₂ 1, MgATP 5, HEPES 10, EGTA 5, pH (KOH) 7.20. After formation of a GΩ seal between the patch electrodes and individual K_v4.3 stably transfected CHO cells, the cell membrane across the pipette tip was ruptured to assure electrical access to the cell interior. All solutions applied to cells were continuously perfused and maintained at room temperature. As soon as a stable seal could be established K_v4.3 currents were measured upon depolarization of the cell membrane from a holding potential of -80 mV to voltages between -120 mV and +60 mV in 20 mV increments of 800 ms duration (Fig. 1). K_v4.3 tail currents were elicited upon a voltage step to 50 mV for 100 ms. The voltage pulses were run at intervals of 5 s.




Fig. 1: Representative K_v4.3 outward currents recorded upon depolarization of the CHO membrane from a holding potential of -80 mV to voltages between -100 mV and +70 mV in 120 mV increments (1000 ms). Voltage pulses to +40 mV (500 ms) elicited K_v4.3 tail currents.

2.2. Biophysical characterization

The K_v4.3 channel is known as an outwardly rectifying potassium channel displaying very fast activation and slow inactivation kinetics. To verify the functional properties of B'SYS's cloned CHO K_v4.3 cells, the IV curve, the activation curve and the inactivation curve were constructed. The voltage dependence of the activation and of the inactivation were best fitted with a Boltzmann function. The half maximal activation was found at -19.9 ± 2.5 mV, $k = 12.7 \pm 0.3$ (n=11), the half maximal inactivation at -49.6 ± 1.4 mV, $k = 5.3 \pm 0.2$ (n=11).




Fig. 2: IV curve of K_v4.3 activation currents. Voltages between -100 and +70 mV were applied.




Fig. 3: Activation (open triangles, blue) and inactivation (closed circles, red) curve of K_v4.3 fitted to Boltzmann function

2.3. KV4.3 channel Blocker Dapoxetine

Dapoxitin, which is known as a K_v4.3 blocker, was used as positive control to validate the pharmacological profile of the CHO K_v4.3 cells. The cells were exposed to 0.53, 5.3 and 53 μ M Dapoxetine and the instantaneous and the steady state current amplitude was analyzed. The inhibition curves (Fig. 5) was best fitted with an IC₅₀ value of 12.4 μ M (Hill Coefficient: 1.2) for the instantaneous current amplitude and an IC₅₀ value of 3.0 μ M (Hill Coefficient: 0.8).




Fig. 4: Representative current recording of a K_v4.3 transfected cell treated with 0.53, 5.3 and 53 μ M Dapoxetine. The cell was depolarized from a holding potential of -80 mV to +40 mV.




Fig. 5: K_v4.3 peak and steady state current inhibition curve of Dapoxetine. The inhibition curve was best fitted with a 3 parameter logistic equation.

2.4. Patch-clamp Success Rates

The patch-clamp properties of the CHO K_v4.3 cell line were elucidated at typical working passage numbers (passage 3 - 15). A total of 27 cells were analyzed. Success for establishment of on-cell configuration was defined as follows: > 1 GΩ. The whole-cell configuration was not accepted if the membrane resistance was below 1 GΩ. A successful recording had to be free of rundown effects and variations in series resistance.

- On-cell successful: **93%** (n=27)
- Whole-cell successful: **89%** (n=27)
- Recording (20-30 min) successful: **78%** (n=27)

3. CELL CULTURE CONDITIONS

3.1. General

CHO K_v4.3 cells are incubated at 37°C in a humidified atmosphere with 5% CO₂ (rel. humidity > 95%). The cells are continuously maintained and passaged in sterile culture flasks containing F12 (HAM) medium supplemented with 10% fetal bovine serum, 1.0% Penicillin/Streptomycin solution and 200 µg/ml Hygromycine. The CHO K_v4.3 cells are passaged at a confluence of about 80%. For electrophysiological measurements the cells are seeded onto e.g. 35 mm sterile culture dishes containing complete medium.

- All solutions and equipment coming in contact with the cells must be sterile.
- Use proper sterile technique and work in a laminar flow hood.
- Be sure to have frozen cell stocks at hand before starting experiments.
- Cells should be split every 2-3 days at 70% - 80% confluence at 1:3 to 1:5 ratio.

3.2. Recommended Complete Medium

- F12 (HAM) with L-Glutamine or GlutaMAX I
- 10% FBS
- 1.0% Penicillin/Streptomycin

3.3. Antibiotics

- CHO K_v4.3 clones were selected under 200 µg/mL Hygromycine antibiotic pressure.
- To cultivate CHO K_v4.3 cells, a reduced antibiotic pressure (200 µg/mL) must be used.
- To separate CHO K_v4.3 cells from untransfected cells, use 500 µg/mL Hygromycine.

Remark: The permanent application of high antibiotic pressure has no effect on current density.

3.4. Thawing Cells

- Remove vial of cells from liquid nitrogen and thaw quickly at 37°C.
- Decontaminate outside of vial with 70% ethanol.
- Transfer cells to a T-75 culture flask containing 10 mL complete medium.
- Incubate cells at 37°C for 4-6 h to allow the cells to attach to the bottom of flask.
- Aspirate off the medium and replace with 10 mL complete medium & antibiotics.
- Antibiotics: 200 µg/mL Hygromycine.
- Incubate cells and check them daily until 70% - 80% confluence is reached.

3.5. Splitting Cells

- When cells are 70% - 80% confluent remove medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipette up and down to break clumps of cells.
- Passage cells into new flask with complete medium and antibiotics at 1:3 to 1:5 ratio.
- Use remaining suspension for counting the cells.

3.6. Freezing Cells

- Cells should have 80% - 90% confluence prior to freezing.
- Remove the complete medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Resuspend Cells in complete medium to reach a density of approximately 1.0E+06 cells per mL.
- Add sterile DMSO to reach 10% DMSO in the complete medium.
- Resuspend cells at a density of approximately 2.0 E+06 cells per mL with fresh freezing medium.
- Mix DMSO and cell suspension by pipetting up and down several times
- Aliquot 0.5 mL of cell suspension into each cryovial.
- Overnight incubate cells in a styropor box at -80°C.
- The next morning transfer cryovial in liquid nitrogen tank for long-term storage.

3.7. Stability of CHO KV4.3 cells

CHO K_v4.3 cells stably express functionally active K_v4.3 potassium channels over 20 passages. Under recommended cell culture conditions no variation in current density was observed during this time.

4. KV4.3 SEQUENCE

4.1. Accession Number NM_002232

Cloned cDNA encodes for the protein of the K_v4.3 channel (NP_004971):

```
MAAGVAAWLPFARAAAIGWMPVANCPMPLAPADKNKRQDELIVLNVSGRRFQTRWTTLER
YDPTLLGSTEKEFFFNEDTKEYFFDRDPEVFRCVLNFYRTGKLHYPRYECISAYDDELAF
YGILPEIIGDCCYEEYKDRKRENAERLMDDNDSENNQESMP SLSFRQTMWRAFENPHTST
LALVFYYVTGFFIAVSVITNVVETVPCGTVPGSKELPCGERYSVAFFCLDTACVMIFTVE
YLLRLFAAPSRYSRFRIRSVMSIIDVVAIMPYYIGLVMTNNEVDVSGAFVTLRVFRVFRIFKF
SRHSQGLRILGYTLKSCASELGFLLFSLTMAIIIFATVMFYAEKGSSASKFTSIPASFWY
TIVTMTTLGYGDMVPKTIAGKIFGSICSLSGVLVIALPVPVIVSNFSRIYHQNRADKRR
AQKKARLARIRVAKTGSSNAYLHSKRNGLLNEALELTGTPEEEHMGKTTSLIESQHHLHL
HCLEKTTGLSYLVDDPLLSVRTSTIKNHEFIDEQMFQNCMESSMQNYPSTRSPSLSSHP
GLTTTCCSRRSCKTTHL PNSLNPATRLRSMQELSTIHIQGSEQPSLTTSRSSLNLKADDG
LRPNCKTSQITTAIISIPTPALTPEGESRPPASP GPNTNIPSIASNVVKVSAL
```

5. CONTACT INFORMATION

5.1. Contact Address for Technical Support and Ordering Information

- B'SYS GmbH
Technology Center Witterswil
Benkenstrasse 254
4108 Witterswil
Switzerland

Tel: +41 61 721 77 44

Fax: +41 61 721 77 41

Email: info@bsys.ch

Web: www.bsys.ch