

B'SYS GmbH

# CHO K<sub>v</sub>4.3 / KChiP2 Cell Line

Specification Sheet

© B'SYS GmbH

## TABLE OF CONTENTS

<b>1. PRODUCT SHIPMENT.....</b>	<b>3</b>
1.1. Product Format.....	3
1.2. Mycoplasma Certificate.....	3
1.3. Assays.....	3
<b>2. VALIDATION OF CHO K<sub>v</sub>4.3 / KChIP2 CELLS .....</b>	<b>4</b>
2.1. Electrophysiology .....	4
2.2. Biophysical characterization .....	5
2.3. K <sub>v</sub> 4.3 channel Blocker Dapoxetine.....	5
2.4. K <sub>v</sub> 4.3 / KChIP2 channel activator NS5806 .....	6
2.5. Automated patch-clamping.....	6
<b>3. CELL CULTURE CONDITIONS.....</b>	<b>7</b>
3.1. General .....	7
3.2. Recommended Complete Medium.....	7
3.3. Antibiotics .....	7
3.4. Thawing Cells .....	7
3.5. Splitting Cells .....	7
3.6. Freezing Cells .....	8
3.7. Stability of CHO K <sub>v</sub> 4.3 /KChIP2 cells .....	8
<b>4. K<sub>v</sub>4.3 AND KChIP2 SEQUENCE.....</b>	<b>9</b>
4.1. Accession Number NM_002232 .....	9
4.2. Accession Number NM_173195 .....	9
<b>5. CONTACT INFORMATION .....</b>	<b>9</b>
5.1. Contact Address for Technical Support and Ordering Information .....	9

## 1. PRODUCT SHIPMENT

### 1.1. Product Format

CHO cells stably transfected with recombinant K<sub>v</sub>4.3 / KChIP2 potassium channel:

- 0.75 mL aliquots of frozen cells at 1 E+06 cells/mL
- Cells are frozen in complete medium with 10% DMSO
- Cells are frozen at passage number, see vial

### 1.2. Mycoplasma Certificate

B'SYS periodically tests cells for presence of mycoplasma by means of highly sensitive PCR based assays. All delivered cells are free of mycoplasma.

### 1.3. Assays

CHO K<sub>v</sub>4.3 KChIP2 cells were validated for manual patch-clamping and automated patch-clamping (Q-Patch, Sophion)

## 2. VALIDATION OF CHO K<sub>v</sub>4.3 / KChIP2 CELLS

### 2.1. Electrophysiology

K<sub>v</sub>4.3 currents were measured by means of the patch-clamp technique in the whole-cell configuration. The bath solution contained (in mM) NaCl 137, KCl 4, CaCl<sub>2</sub> 1.8, MgCl<sub>2</sub> 1, D-glucose 10, HEPES 10, pH (NaOH) 7.40. The pipette solution consisted of (in mM) KCl 130, MgCl<sub>2</sub> 1, MgATP 5, HEPES 10, EGTA 5, pH (KOH) 7.20. After formation of a GΩ seal between the patch electrodes and individual K<sub>v</sub>4.3 / KChIP2 stably transfected CHO cells, the cell membrane across the pipette tip was ruptured to assure electrical access to the cell interior. All solutions applied to cells were continuously perfused and maintained at room temperature. As soon as a stable seal could be established K<sub>v</sub>4.3 / KChIP2 currents were measured upon depolarization of the cell membrane from a holding potential of -80 mV to voltages between -120 mV and +60 mV in 20 mV increments of 800 ms duration (Fig. 1). K<sub>v</sub>4.3 / KChIP2 tail currents were elicited upon a voltage step to 50 mV for 100 ms. The voltage pulses were run at intervals of 5 s.





Fig. 1: Representative K<sub>v</sub>4.3 (black) and K<sub>v</sub>4.3 / KChIP2 (red) outward currents recorded upon depolarization of the CHO membrane from a holding potential of -80 mV to voltages between -100 mV and +70 mV in 20 mV increments (1000 ms). Voltage pulses to +40 mV (500 ms) elicited K<sub>v</sub>4.3, K<sub>v</sub>4.3 / KChIP2 tail currents.


## 2.2. Biophysical characterization

To differentiate between K<sub>v</sub>4.3 and K<sub>v</sub>4.3 / KChIP2 the kinetic of inactivation was analyzed by fitting a monoexponential function to the inactivated current observed during depolarizations to +60 mV. The time constant for channels co-expressed with KChIP2 is about twice slower compared to the K<sub>v</sub>4.3 channel alone.




## 2.3. K<sub>v</sub>4.3 channel Blocker Dapoxetine

Dapoxetine, which is known as a K<sub>v</sub>4.3 blocker, was used as positive control to validate the pharmacological profile of the CHO K<sub>v</sub>4.3 cells. The cells were exposed to 0.53, 5.3 and 53 μM Dapoxetine and the instantaneous and the steady state current amplitude was analyzed. The inhibition curves (Fig. 3) were best fitted with an IC<sub>50</sub> value of 12.4 μM (Hill Coefficient: 1.2) for the instantaneous current amplitude and an IC<sub>50</sub> value of 3.0 μM (Hill Coefficient: 0.8).




**Fig. 2:** Representative current recording of a K<sub>v</sub>4.3 transfected cell treated with 0.53, 5.3 and 53 μM Dapoxetine. The cell was depolarized from a holding potential of -80 mV to +40 mV.



**Fig. 3:** K<sub>v</sub>4.3 peak and steady state current inhibition curve of Dapoxetine. The inhibition curve was best fitted with a 3 parameter logistic equation.

## 2.4. K<sub>v</sub>4.3 / KChIP2 channel activator NS5806


NS5806 (1-[2,4-dibromo-6-(1H-tetrazol-5-yl)-phenyl]-3-(3,5-bis-trifluoromethyl-phenyl)-urea), known as a K<sub>v</sub>4.3 activator, was also as positive control to validate the pharmacological profile of the CHO K<sub>v</sub>4.3 / KChIP2 cells. The cells were exposed to 10 μM NS5806 and the instantaneous and the steady state current amplitudes were analyzed. The application of 10 μM NS5806 induced an increase of K<sub>v</sub>4.3 / KChIP2 peak current amplitudes by about 10% and slowed down the inactivation significantly (P=0.048).



**Fig. 4:** Representative current recording of a K<sub>v</sub>4.3 transfected cell treated with 10 μM NS5806. The cell was depolarized from a holding potential of -80 mV to +40 mV. Black: Control, Red: 10 μM NS5806

## 2.5. Automated patch-clamping

CHO K<sub>v</sub>4.3 / KChIP2 cells were tested on the Q-Patch HTX and Q-Patch 16X. The cells showed current amplitudes larger than 10 nA when depolarized from -80 mV to +40 mV.



**Fig. 5:** Automated patch clamping: left: representative current recording of a K<sub>v</sub>4.3 / KChIP2 transfected cell treated with increasing concentrations of Quinidine. The cell was depolarized from a holding potential of -80 mV to +40 mV. right: Statistic of success of an average Q-Plate tested on the Q-Patch16X (one column).

## 3. CELL CULTURE CONDITIONS

### 3.1. General

CHO K<sub>v</sub>4.3 / KChIP2 cells are incubated at 37°C in a humidified atmosphere with 5% CO<sub>2</sub> (rel. humidity > 95%). The cells are continuously maintained and passaged in sterile culture flasks containing F12 (HAM) medium supplemented with 10% fetal bovine serum, 1.0% Penicillin/Streptomycin solution and 200 µg/mL Hygromycine and 10 µg/mL Puromycin. The CHO K<sub>v</sub>4.3 / KChIP2 cells are passaged at a confluence of about 80%. For manual electrophysiological measurements the cells are seeded onto e.g. 35 mm sterile culture dishes containing complete medium.

- All solutions and equipment coming in contact with the cells must be sterile.
- Use proper sterile technique and work in a laminar flow hood.
- Be sure to have frozen cell stocks at hand before starting experiments.
- Cells should be split every 2-3 days at 70% - 80% confluence at 1:3 to 1:5 ratio.

### 3.2. Recommended Complete Medium

- F12 (HAM) with L-Glutamine or GlutaMAX I
- 10% FBS
- 1.0% Penicillin/Streptomycin

### 3.3. Antibiotics

- CHO K<sub>v</sub>4.3 / KChIP2 clones were selected under 200 µg/mL Hygromycine antibiotic pressure.
- To cultivate CHO K<sub>v</sub>4.3 / KChIP2 cells, a reduced antibiotic pressure (200 µg/mL Hygromycine and 10 µg/mL Puromycin) must be used.
- To separate CHO K<sub>v</sub>4.3 / KChIP2 cells from untransfected cells, use 500 µg/mL Hygromycine and 15 µg/mL Puromycin.

Remark: The permanent application of high antibiotic pressure has no effect on current density.

### 3.4. Thawing Cells

- Remove vial of cells from liquid nitrogen and thaw quickly at 37°C.
- Decontaminate outside of vial with 70% ethanol.
- Transfer cells to a T-75 culture flask containing 10 mL complete medium.
- Incubate cells at 37°C for 4-6 h to allow the cells to attach to the bottom of flask.
- Aspirate off the medium and replace with 10 mL complete medium & antibiotics.
- Antibiotics: 200 µg/mL Hygromycine and 10 µg/mL Puromycin.
- Incubate cells and check them daily until 70% - 80% confluence is reached.

### 3.5. Splitting Cells

- When cells are 70% - 80% confluent remove medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipette up and down to break clumps of cells.
- Passage cells into new flask with complete medium and antibiotics at 1:3 to 1:5 ratio.
- Use remaining suspension for counting the cells.

### 3.6. Freezing Cells

- Cells should have 80% - 90% confluence prior to freezing.
- Remove the complete medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Resuspend cells in complete medium to reach a density of approximately 1.0E+06 cells per mL.
- Add sterile DMSO to reach 10% DMSO in the complete medium.
- Resuspend cells at a density of approximately 2.0 E+06 cells per mL with fresh freezing medium.
- Mix DMSO and cell suspension by pipetting up and down several times
- Aliquot 0.5 mL of cell suspension into each cryovial.
- Overnight incubate cells in a styropor box at -80°C.
- The next morning transfer cryovial in liquid nitrogen tank for long-term storage.

### 3.7. Stability of CHO K<sub>v</sub>4.3 /KChIP2 cells

CHO K<sub>v</sub>4.3 cells stably express functionally active K<sub>v</sub>4.3 / KChIP2 potassium channels over 20 passages. Under recommended cell culture conditions no variation in current density was observed during this time.



## 4. K<sub>v</sub>4.3 AND KCHIP2 SEQUENCE

### 4.1. Accession Number NM\_002232

Cloned cDNA encodes for the protein of the long (cardial) version of the K<sub>v</sub>4.3 channel (NP\_004971) protein:

```
MAAGVAAWLPFARAAAIGWMPVANCPMLAPADKNKRQDELIVLNVSGRRFQTRWTTLER
YPTLLGSTEKEFFFNEDTKEYFFDRDPEVFRVCLNLFYRTGKLLHYPRYECISAYDDELAF
YGILPEIIGDCCYEEYKDRKRENAERLMDNDSENNQESMPSLSFRTMWRAFENPHTST
LALVFYVYTGFFIAVSVITNVVETVPCGTVPGSKELPCGERYSVAFFCLDTACVMIFTVE
YLLRFLFAAPSRYRFIRSVMSIIDVVAIMPYYIGLVMTNNEVDVSGAFVTLRVFRVFRIFKF
SRHSQGLRILGYTLKSCASELGFLLFSLTMAIIIFATVMFYAEKGSSASKFTSIPASFWY
TIVTMTTLGYGDMVPKTIAGKIFGSI CSLSGVLVIALPVPVIVSNFSRIYHQNRADKRR
AQKKARLARIRVAKTGSSNAYLHSCRNGLLNEALELTGTPEEEHMGKTTSLIESQHHLHL
HCLEKTTGLSYLVDDPLLSVRTSTIKNHEFIDEQMFEQNCMESSMQNYPSTRSPSLSSHP
GLTTTCCSRRSKKTTHLPNSNLPATRLRSMQELSTIHIQGEQPSLTTSRSSLNLKADDG
LRPNCKTSQITTAIISIPTPPALTPEGESRPPASPGPNTNIPSIASNVVKVSAL
```

### 4.2. Accession Number NM\_173195

Cloned cDNA encodes for the protein of the KChIP2 (NP\_775287) protein:

```
MRGQGRKESLSDSRDLGSDYDQLTDSVDDEFELSTVCHRPEGLEQLQEQTFRKELQVL
YRGFKNECPSGI VNEENFKQIYSQFFPQGDSSYATFLFNAFDTNHDGVSFEDFVAGLS
VILRGTVDDRNLNWFNLYDLNKDGCITKEEMLDIMKSIYDMMGKYTYPALREEAPREHVE
SFFQKMDRNDKGVVTIEEFIESQKDENIMRSMQLFDNVI
```

## 5. CONTACT INFORMATION

### 5.1. Contact Address for Technical Support and Ordering Information

- B'SYS GmbH  
Technology Center Witterswil  
Benkenstrasse 254  
4108 Witterswil  
Switzerland

Tel: +41 61 721 77 44  
Fax: +41 61 721 77 41  
Email: [info@bsys.ch](mailto:info@bsys.ch)  
Web: [www.bsys.ch](http://www.bsys.ch)