

B'SYS GmbH

# CHO K<sub>v</sub>7.2/3 Cell Line

Specification Sheet

© B'SYS GmbH

---

## TABLE OF CONTENTS

<b>1. BACKGROUND</b> .....	<b>3</b>
1.1. The neuronal M-current .....	3
1.2. B'SYS CHO K <sub>v</sub> 7.2 and CHO K <sub>v</sub> 7.2/7.3 Cells .....	3
<b>2. PRODUCT SHIPMENT</b> .....	<b>4</b>
2.1. Product Format .....	4
2.2. Mycoplasma Certificate.....	4
<b>3. VALIDATION OF CHO K<sub>v</sub>7.2 AND K<sub>v</sub>7.2/7.3</b> .....	<b>5</b>
3.1. Electrophysiology .....	5
3.2. Positive Control XE-991 .....	6
<b>4. HETEROMERIC K<sub>v</sub>7.2/7.3</b> .....	<b>7</b>
4.1. Current amplitudes of Homo- and Heteromers .....	7
4.2. Positive control: TEA .....	7
4.3. Potassium channel activator: Retigabine .....	8
<b>5. CELL CULTURE CONDITIONS</b> .....	<b>9</b>
5.1. General.....	9
5.2. Recommended Complete Medium .....	9
5.3. Antibiotics .....	9
5.4. Thawing Cells .....	9
5.5. Splitting Cells .....	10
5.6. Freezing Medium.....	10
5.7. Freezing Cells .....	10
5.8. Stability of CHO K <sub>v</sub> 7.2 and K <sub>v</sub> 7.2/7.3 cells.....	10
<b>6. SEQUENCES</b> .....	<b>11</b>
6.1. Human K <sub>v</sub> 7.2 Accession Number NM172107.2 .....	11
6.2. Human K <sub>v</sub> 7.3 Accession Number NM001204824 .....	11
<b>7. CONTACT INFORMATION</b> .....	<b>11</b>

---

## 1. BACKGROUND

### 1.1. The neuronal M-current

K<sub>v</sub>7.2 and K<sub>v</sub>7.3 are voltage-gated potassium (K<sup>+</sup>) channels, belonging to the subfamily of K<sub>v</sub>7.x. They are encoded by *KCNQ2* and *KCNQ3* genes respectively; and co-assemble together to form an M-type homomeric or heteromeric K<sub>v</sub> channel. These channels, having the properties of the native M-current, is a slowly activating and deactivating (delayed rectifier) K<sup>+</sup> channel that regulates synaptic plasticity and neuronal excitability by causing delayed membrane hyperpolarization.

K<sub>v</sub>7.2 and 7.3 are found in hippocampus, cortex, thalamus, cerebellum and brain stem. Further they are found in nodes of Ranvier, sympathetic and dorsal root ganglia. Loss of function mutations in *KCNQ* genes can cause repeated neuronal discharges and is associated with rare hereditary human epilepsy called benign familial neonatal convulsions (BFNC). Furthermore, *KCNQ* genes are expressed in sensory neurons and are therefore also involved in nociceptive signaling.

K<sub>v</sub>7.2/7.3 inhibitors are developed for the treatment of learning and memory disorders; whereas K<sub>v</sub>7.2/7.3 activators for the treatment of epilepsy and pain.

### 1.2. B'SYS CHO K<sub>v</sub>7.2 and CHO K<sub>v</sub>7.2/7.3 Cells

As a new test system B'SYS has designed a CHO K<sub>v</sub>7.2 cell line with constitutive expression of K<sub>v</sub>7.2 channels. The K<sub>v</sub>7.2 cDNA was cloned and transfected into CHO cells and then the functional properties of the K<sub>v</sub>7.2 channels validated by means of the patch-clamp technique.

The K<sub>v</sub>7.2 and K<sub>v</sub>7.3 cDNAs were cloned in the same bidirectional expression vector and were transfected on CHO cells. Their functional properties are also validated by the means of the manual patch-clamp technique.

## 2. PRODUCT SHIPMENT

### 2.1. Product Format

CHO cells stably transfected with recombinant human K<sub>v</sub>7.2 or K<sub>v</sub>7.2/7.3 ion channel:

- 1 x 0.8 ml aliquots of frozen cells at 2.3 E+06 cells/mL
- Cells are frozen in complete medium with 10 % DMSO


### 2.2. Mycoplasma Certificate

B'SYS periodically tests cells for presence of mycoplasma by means of highly sensitive PCR based assays. All delivered cells are free of mycoplasma.

### 3. VALIDATION OF CHO K<sub>V</sub>7.2 AND K<sub>V</sub>7.2/7.3


#### 3.1. Electrophysiology

K<sub>V</sub>7.2 and K<sub>V</sub>7.2/7.3 currents were measured by means of the patch-clamp technique in the perforated-patch configuration. Access to cells was established with Amphotericin B which perforates the patched membrane. The bath solution contained (in mM) Potassium Chloride 137, Potassium Chloride 4, Calcium Chloride 1.8, Magnesium Chloride 1, HEPES 10, D-Glucose 10, pH (NaOH) 7.4. The pipette solution consisted of (in mM) Potassium Chloride 130, Magnesium Chloride 1, Mg-ATP 5, HEPES 10, EGTA 5, Amphotericin B (0.5 mg/mL), pH (KOH) 7.2. After formation of a GΩ seal between the patch electrodes and individual K<sub>V</sub>7.2 or K<sub>V</sub>7.2/7.3 stably transfected CHO cells, the cell membrane was perforated within 5 minutes. Improvement of electrical access to the cell interior was monitored until series resistance was stable. All solutions applied to cells were continuously perfused and maintained at room temperature.



### 3.2. Positive Control XE-991


XE-991, which is known a K<sub>V</sub> blocker, was used as positive control to validate the pharmacological profile of the CHO cells. The effect of XE-991 was first measured at 10 μM, which blocked K<sub>V</sub>7.2 currents almost completely ( $10.97 \pm 3.06\%$  relative current amplitude, mean  $\pm$  SEM of 5 cells). When the cells were exposed to 0.3, 1.0, 3.0 and 10 μM XE-991, the steady-state level relative current amplitudes revealed values of  $89.46 \pm 3.88\%$ ,  $52.93 \pm 7.41\%$ ,  $27.53 \pm 5.96$  and  $10.97 \pm 3.06\%$ , respectively (mean  $\pm$  SEM of 5 cells). The inhibition curve (Fig. 2B) was best fitted with an IC<sub>50</sub> value of **1.27 μM**. Similarly, for the heteromeric K<sub>V</sub>7.2/7.3 currents were blocked by XE-991. When the cells were exposed to 0.1, 1.0, 3.0 and 10 μM XE-991, the steady-state level relative current amplitudes revealed values of  $83.64 \pm 6.43\%$ ,  $61.74 \pm 4.21\%$ ,  $28.38 \pm 4.87$  and  $15.62 \pm 2.41\%$ , respectively (mean  $\pm$  SEM of 5 cells). The inhibition curve (Fig. 2D) was best fitted with an IC<sub>50</sub> value of **1.29 μM**.



## 4. HETEROMERIC Kv7.2/7.3


### 4.1. Current amplitudes of Homo- and Heteromers

For Kv7.2/7.3 heteromeric channels, the current amplitude was strongly reduced in comparison to Kv7.2 channels, but still larger than the current amplitude that was observed in cells that were expressing with Kv7.3 only.



### 4.2. Positive control: TEA

Kv7.2 currents could be blocked by 10  $\mu$ M TEA. The remaining current amplitude was 8%. Cells transfected with Kv7.2 / 7.3 were less sensitive, only 65% of the current amplitude was blocked.



### 4.3. Potassium channel activator: Retigabine

Retigabine is a novel anti-convulsant known to act as a positive modulator of the heteromer K<sub>v</sub>7.2/7.3. At a concentration of 10  $\mu$ M the current for the K<sub>v</sub>7.2/7.3 was stimulated by 63%, while the two homomers were only stimulated by 30%. For K<sub>v</sub>7.2/7.3, the dose response curve was generated and the EC<sub>50</sub> value was calculated to be 8.34  $\mu$ M.




Fig.5 Effect of Retigabine on K<sub>v</sub>7.2/7.3 channels. Top left: Representative current recording of a K<sub>v</sub>7.2/7.3 transfected cell perfused with increasing concentrations of Retigabine. Top right: dose response curve of Retigabine for K<sub>v</sub>7.2/7.3 channels. Bottom: Summary of current stimulation with 10  $\mu$ M Retigabine.



---

## 5. CELL CULTURE CONDITIONS

### 5.1. General

CHO K<sub>v</sub>7.2 and K<sub>v</sub>7.2/7.3 cells are incubated at 37°C in a humidified atmosphere with 5% CO<sub>2</sub> (rel. humidity > 95%). The cells are continuously maintained and passaged in sterile culture flasks containing F12 (HAM) medium supplemented with 10% fetal bovine serum, 1% Penicillin/Streptomycin solution and 500 µg/mL Hygromycin (K<sub>v</sub>7.2) or 5 µg/mL Puromycin (K<sub>v</sub>7.2/7.3) and 100 µg/mL Zeocin. The CHO cells are passaged at a confluence of about 50-80%. For electrophysiological measurements the cells are seeded onto e.g. 35 mm sterile culture dishes containing complete medium.

- All solutions and equipment coming in contact with the cells must be sterile.
- Use proper sterile technique and work in a laminar flow hood.
- Be sure to have frozen cell stocks at hand before starting experiments.
- Cells should be split every 2-3 days at 70% - 80% confluency at 1:3 to 1:5 ratio.

### 5.2. Recommended Complete Medium

- F12 (HAM) with L-Glutamine or GlutaMAX I
- 10% FBS
- 1% Penicillin/Streptomycin

### 5.3. Antibiotics

- CHO K<sub>v</sub>7.2 clones were selected under 1000 µg/mL Hygromycin antibiotic pressure.
- To cultivate CHO K<sub>v</sub>7.2 cells, also a reduced antibiotic pressure (500 µg/mL) can be used.
- To separate CHO K<sub>v</sub>7.2 cells from untransfected cells, use 1000 µg/mL Hygromycin.
  
- CHO K<sub>v</sub>7.2/7.3 clones were selected under 10 µg/mL Puromycin and 100 µg/mL Zeocin antibiotic pressure.
- To cultivate CHO K<sub>v</sub>7.2/7.3 cells, also a reduced antibiotic pressure (5 µg/mL Puromycin and 100 µg/mL Zeocin) can be used.
- To separate CHO K<sub>v</sub>7.2/7.3 cells from untransfected cells, use 10 µg/mL Puromycin and 100 µg/mL Zeocine
- 

Remark: The permanent application of high antibiotic pressure has no effect on current density.

### 5.4. Thawing Cells

- Remove vial of cells from liquid nitrogen and thaw quickly at 37°C.
- Decontaminate outside of vial with 70% ethanol.
- Transfer cells to a T-75 culture flask containing 10 mL complete medium.
- Incubate cells at 37°C for 4-6 h to allow the cells to attach to the bottom of flask.
- Aspirate off the medium and replace with 10 mL complete medium & antibiotics.
- Antibiotics: 500 µg/mL Hygromycin (K<sub>v</sub>7.2); 5 µg/mL Puromycin and 100 µg/mL Zeocine (K<sub>v</sub>7.2/7.3).
- Incubate cells and check them daily until 70% - 80% confluency is reached.

### 5.5. Splitting Cells

- When cells are 70% - 80% confluent remove medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Passage cells into new flask with complete medium and antibiotics at 1:3 to 1:5 ratio.
- Use remaining suspension for counting the cells.

### 5.6. Freezing Medium

- Mix 0.9 mL fresh complete medium and 0.1 mL DMSO for every 1 mL freezing medium.
- Sterilize freezing medium by means of appropriate micro filter (0.1 µm – 0.2 µm).

### 5.7. Freezing Cells

- Prepare fresh freezing medium and keep it on ice.
- Cells should have 80% - 90% confluency prior to freezing.
- Remove the complete medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Pellet cells with centrifuge and carefully aspirate off medium.
- Resuspend cells at a density of approximately 2.0 E+06 cells per mL with fresh freezing medium.
- Aliquot 0.8 mL of cell suspension into each cryovial.
- Overnight incubate cells in a styropor box at –80°C.
- The next morning transfer cryovial in liquid nitrogen tank for long-term storage.

### 5.8. Stability of CHO Kv7.2 and Kv7.2/7.3 cells

CHO Kv7.2 and Kv7.2/7.3 cells stably express functionally active Kv7.2 or Kv7.2/7.3 potassium channels over 20 passages. Under recommended cell culture conditions no variation in current density was observed over 20 cell splitting cycles.

## 6. SEQUENCES

Human K<sub>v</sub>7.2 Accession Number NP\_742105.1  
Human K<sub>v</sub>7.3 Accession Number NP\_004510.1

## 7. CONTACT INFORMATION

- B'SYS GmbH  
Technology Center Witterswil  
Benkenstrasse 254  
4108 Witterswil  
Switzerland

Tel: +41 61 721 77 44  
Fax: +41 61 721 77 41  
Email: [info@bsys.ch](mailto:info@bsys.ch)  
Web: [www.bsys.ch](http://www.bsys.ch)